



Division of
**Consolidated
Laboratory Services**



DEPARTMENT OF
GENERAL SERVICES



Virginia Bioscience Foundation in partnership with
Varina High School and the Division of
Consolidated Laboratory Services

Bacteria Laboratory Presentation

January 23, 2024







1

Division of
**Consolidated
Laboratory
Services**









DEPARTMENT OF
GENERAL SERVICES

Tentative Schedule

Time	Activity
10:00 a.m.	DCLS Arrival and Set Up
10:50 a.m.	Introductions
11:00 a.m.	Station 1: Overview of the Lab (All Students)
11:20 a.m.	Begin 1 st Station Assignment (2,3,4 or 5)
11:40 a.m.	Begin 2 nd Station (2,3,4 or 5)
12:00 p.m.	Begin 3 rd Station (2,3,4, or 5)
12:25 p.m.	Scheduled HCPS Break
12:30 p.m.	Begin Last Station Assignment (2,3,4, or 5)
12:50 p.m.	Close-out and Clean-up

2



Special Thanks!

- Virginia Bioscience Foundation
 - Eric Rhoades
 - Jim Powers
- Go Virginia Grant
- Varina Center for Environmental Studies and Sustainability
 - Jill Rich
 - Kassandra Epps
 - Cherita Sears
- Varina Center for Communications and Media
- Department of General Services
 - Division of Consolidated Laboratory Services
 - Office of Communications
- Collaborative Teaching and Learning Group, LLC

3

4



5

Definitions

- Ampule: a sealed glass container containing a liquid or powder
- Aseptic Technique: Procedures used to prevent contamination
- Binary Fission: a process where an organism duplicates its genetic material and then divides into two parts which each new organism “part” receiving one copy of the genetic material
- Biological Metabolism: a series of reactions that occur within living organisms to sustain life
- Differentiation: the act of showing or finding difference between things that are compared (distinguish one microbe from another)
- Don: To put on
- Doff: To take off or remove



6

Definitions

- Incubate: to maintain an item for a defined period of time under known or well-defined conditions (ie. temperature)
- Infectious: The ability to cause disease
- Impregnated: Soaked or saturated with a substance
- Lawn of Bacteria: a continuous, uniform layer of bacterial growth on the surface of agar media
- Emulsify: to mix so thoroughly that it becomes an emulsion, or a mixture of two substances that can't be completely blended together
- Homogenous: of uniform or the same structure or composition throughout



Definitions

- McFarland standards: a series of standards of varying turbidities that allow scientist to compare unknowns to visually to approximate the number of bacteria in a liquid suspension
- Media: liquid or gel that contains nutrients to help organisms grow
- Microbe: a very small living organism that cannot be seen without a microscope
- Morphology: the form and structure of living organisms such as bacteria
- Phenotypic: an observable trait

7

**Bacteria and Virus
Station 1: Overview**

Presenters:
January 23, 2024

8



9

Learning Objectives

- Practice effective safety protocols during all stations
- Discuss principles of bacteriology and virology
- Examine laboratory techniques for the identification and differentiation of microbes



10

Procedural Steps

- Meet your facilitators
- Divide into groups
- Rotate to each station
- Follow station facilitator instructions
- Complete the hands-on activity
- Ask questions
- Record observations and reflections
- Provide feedback to facilitators, teachers, or other guests



11

Safety

- **Universal precautions:** handle all materials as potentially infectious
- Know the location of eye wash/drench hose
- Wear fully enclosed shoes (heel and toe)
- **Personal Protective Equipment (PPE):** Don (which means to put on) required PPE including, lab coat, safety glasses, and gloves
- Understand hazards associated with chemicals and biological organisms described at each station
- **No food, drink or gum chewing** at stations
- **Wash hands** thoroughly with soap and water



12

Overview of Stations

- Station 1: Overview
- Station 2: Culture and Sensitivity
- Station 3: Gram Stain
- Station 4: Analytical Profile Index (API) Set-Up
- Station 5: Analytical Profile Index (API) Interpretation



13

Materials

- Station instructions
- Videos
- Writing utensil (permanent ink not pencil)
- Personal protective equipment (PPE)
 - Lab coat
 - Safety glasses
 - Nitrile gloves



14

Procedural Steps

- Station 2: Culture and Sensitivity
 - Bacterial cultures
 - Growth media
 - Aseptic technique
 - Colony morphology
 - Size, texture, transparency, pigmentation, elevation, form
 - Antimicrobial resistance



Procedural Steps

- Station 3: Gram Stain
 - Cell walls
 - Chemical and physical properties
 - Gram stain for differentiation
 - Microscopic examination
 - Antibiotics

15



Procedural Steps

- Station 4: Analytical Profile Index (API) Set-Up
 - Aseptic technique
 - Pipetting
- Station 5: Analytical Profile Index (API) Interpretation
 - Phenotypic microbial identification
 - Biochemical metabolism
 - Color changes
 - Metabolic byproducts (i.e. gases, enzymatic activity)

16



17

Resources

- Antimicrobial resistance
 - <https://www.cdc.gov/drugresistance/index.html>
- Colony characteristics
 - https://www.cdc.gov/labtraining/docs/job_aids/biochemicals_gram_positive_organism_id/Colonial_Characteristics_Branded_508.pdf
- Gram stain
 - https://www.cdc.gov/labtraining/docs/job_aids/biochemicals_and_gram-negative_organism_id/Gram_Stain_Procedure_Branded_508.pdf



18

Resources

- Culture media
 - https://www.cdc.gov/labtraining/docs/job_aids/biochemicals_and_gram-negative_organism_id/Culture_Media_Table_Branded_508.pdf
- Microbiological identification
 - <https://www.cdc.gov/microbenet/about.html>
- API
 - <https://phil.cdc.gov/Details.aspx?pid=19231>

DCLS

Division of
**Consolidated
Laboratory Services**

DGS
DEPARTMENT OF
GENERAL SERVICES

Bacteriology and Virology

Station 2:
Culture and Sensitivity

Presenters:

January 23, 2024

19

Division of
**Consolidated
Laboratory
Services**

DGS
DEPARTMENT OF
GENERAL SERVICES

Learning Objectives

- Discuss phenotypic characteristics of bacteria
- Prepare cultures using baker's yeast, not bacteria
- Observe previously prepared bacterial plates to compare phenotypic characteristics and record observations
- Recognize resistance patterns for sensitive versus resistant bacteria

20




Bacterial Cultures

- Simple cell division by binary fission
- Bacteria multiply once every 20-30 minutes
- Logarithmic growth – in 1 hour a single bacterium can produce 8 or more bacteria
- Aseptic technique
 - Procedures used to prevent contamination
- Growth media (agar)
 - Contents: carbohydrates, minerals, vitamins, proteins

21



Binary Fission



22



23

Colony Morphology

- Texture – appearance of the surface
 - Smooth, shiny, mucoid, dry, flaky, etc.
- Transparency – see through them
 - Transparent, translucent, opaque
- Color or pigmentation – intracellular pigments
 - Yellow, pink, purple, red, etc.


























24

Colony Morphology

- Size of Colony
 - Punctiform, small, medium, large
- Elevations
 - Flat, raised, convex, pulvinate, umbonate
- Forms
 - Circular, rhizoid, irregular, filamentous
- Hemolysis
 - Beta, alpha, gamma



Colony Morphology

shape	size	surface	color	opacity	elevation	margin
 circular	 small	smooth	 white	transparent	 flat	 even
 punctiform	 medium	glistening	 creamy-white	translucent	 umbonate	 wavy
 filamentous	 large	rough	 yellow	opaque	 raised	 filamentous
 irregular		wrinkle	 orange		 convex	 lobate
 rhizoid		dull	 green		 pulvinate	 curled

rsscience.com

25




Safety

- **Universal precautions:** handle all materials as potentially infectious
- Know the location of eye wash/drench hose
- Wear fully enclosed shoes (heel and toe)
- **Personal Protective Equipment (PPE):** Don (which means to put on) required PPE including fully enclosed shoes, lab coat, safety glasses, and gloves
- Understand hazards associated with chemicals and biological organisms described at each station
- **No food, drink or gum chewing** at stations
- **Wash hands** thoroughly with soap and water

26



Inoculation Demonstration



27



Materials and Reagents

- Gloves
- Inoculation loops
- McFarland standards
- Blood agar plates
- MacConkey plates
- Mueller Hinton plates
- Baker's yeast

28



29

Procedural Steps

- Don appropriate PPE
- Obtain yeast sample
- Using aseptic technique, retrieve some of the yeast sample with the disposable inoculating loop
- Place the loop flat on the agar surface at the top of the plate and drag in a zig zag pattern for approximately $\frac{1}{4}$ of the plate
- Discard the used loop and retrieve a new loop





30

Procedural Steps

- Place the new loop flat on the agar surface in the quadrant that was previously streaked
- Drag the loop in a zigzag pattern for $\frac{1}{4}$ of the plate
- Discard loop and retrieve a new one
- Repeat this process until all 4 quadrants have been streaked



Four Quadrant Streak Plate Method



31



Colony Description: MacConkey Agar

Escherichia coli

- Size: _____
- Texture: _____
- Transparency: _____
- Pigmentation: _____
- Elevation: _____
- Form: _____
- Hemolysis: _____

Pseudomonas aeruginosa

- Size: _____
- Texture: _____
- Transparency: _____
- Pigmentation: _____
- Elevation: _____
- Form: _____
- Hemolysis: _____

32



33

Colony Description: Blood Agar

- | | |
|--|---|
| <ul style="list-style-type: none"> • <i>Staphylococcus aureus</i> • Size: _____ • Texture: _____ • Transparency: _____ • Pigmentation: _____ • Elevation: _____ • Form: _____ • Hemolysis: _____ | <ul style="list-style-type: none"> • <i>Escherichia coli</i> • Size: _____ • Texture: _____ • Transparency: _____ • Pigmentation: _____ • Elevation: _____ • Form: _____ • Hemolysis: _____ |
|--|---|



34

Antimicrobial Resistance Kirby Bauer Disk Diffusion

- Determines whether an organism is sensitive or resistant to an antibiotic
- Used by physicians to make treatment decisions
- Need a "lawn" of bacterial growth and discs impregnated with antibiotics
- Size of the zone of inhibition determines whether a bacteria is sensitive or resistant to that antibiotic
- *Note: different test methods are used for environmental testing*




Antimicrobial Resistance Kirby Bauer Disk Diffusion

Sensitive



Resistant



35




Antimicrobial Resistance -Kirby Bauer Disk Diffusion



36



Antimicrobial Resistance Kirby Bauer Disk Diffusion



37



Division of
**Consolidated
Laboratory Services**



DEPARTMENT OF
GENERAL SERVICES




Bacteriology and Virology

Station 3:
Gram Stain

Presenters:
January 23, 2024

38



39

Learning Objective

- Describe the reason for performing gram stains
- Explain 3 outcomes from staining cells and whether they categorize as gram positive, negative, or variable
- Understand the reason why cells stain purple/blue vs red/pink
- Summarize steps for completing a gram stain
- Explain modes of action for antibiotics



40


Purpose

- A differential staining procedure to characterize the chemical and physical properties of a cell wall
- Used to differentiate bacteria based on cellular form, size, morphology, and gram reaction
- Guides the laboratory in the selection of appropriate isolation media and tests needed for identification
- Critical test and in some case can serve as the presumptive diagnosis for infectious agents



Bacterial Characteristics


- Size
 - 1 to 10 microns – microscopic range
- Morphology



41



Bacterial Cell Wall Characteristics




42



Antibiotics

- Medicines that can kill bacteria *OR* prevent them from multiplying
- 1st Natural Antibiotic – Penicillin (1928)
- Different antibiotics have different modes of action to disrupt bacterial cell growth
 - Interfere with reproduction
 - Block protein production
 - Attack cell wall




43



Modes of Resistance


- Many environmental bacteria are multi-drug resistant
- Target modification - Evolution and genetic mutation
- Efflux - Pump out the antibiotics
- Bypass and enzyme-catalyzed destruction - genetic, enzymatic, or replacement of the original target (bypass)



44



Antibiotics



45



Safety

- **Universal precautions:** handle all materials as potentially infectious
- Know the location of eye wash/drench hose
- Wear fully enclosed shoes (heel and toe)
- **Personal Protective Equipment (PPE):** Don (which means to put on) required PPE including fully enclosed shoes, lab coat, safety glasses, and gloves
- Understand hazards associated with chemicals and biological organisms described at each station
- **No food, drink or gum chewing** at stations
- **Wash hands** thoroughly with soap and water

46



Materials

- Gloves
- Lens paper
- Light microscope
- Prepared gram stain slides
- Immersion oil
- Pasteur pipettes
- Alcohol swabs




47



Procedural Steps

1. Prepare a thin smear of the material for study on a microscope slide

- Place 1-2 drops of broth (organism suspended in media) on a slide using a sterile pipette
- Use a slide warmer to heat fix the smear




48



Procedural Steps

2. Stain the smear

- Place the heat fixed smear on a staining rack over the sink and overlay with Gram Crystal Violet for 1 minute
- Rinse with water and overlay with reconstituted Gram Iodine mordant for 1 minute
- Flood with Gram Decolorizer until the solvent flows colorless from the slide (10-30 sec)




49



Procedural Steps

- Rinse with water and overlay with Gram Safranin for 30 seconds
- Rinse with water and allow to dry
- Examine the slide microscopically with the 100X oil immersion objective



50





Procedural Steps - Video




51




Procedural Steps

3. Interpretation

Gram positive organisms stain deep blue or purple





Gram negative organisms stain red or pink.



Gram variable organisms stain a combination of purple and pink



52


 DCLS
 Division of
**Consolidated
 Laboratory Services**

 DEPARTMENT OF
 GENERAL SERVICES

Bacteriology and Virology Station 4: Analytical Profile Index Set-Up

Presenters:
 January 23, 2024

53

Division of
**Consolidated
 Laboratory
 Services**


 DEPARTMENT OF
 GENERAL SERVICES

Learning Objectives

- Define the purpose of an Analytical Profile Index (API) strip
- Recall different types of organisms that can be identified using API
- Describe biochemicals
- Summarize the steps for preparing an API strip

54



55

Purpose

Rapid and accurate species identification of known, clinically-relevant microorganisms

Gram negative organisms:

- *Neisseria* sp.
- *Haemophilus*
- *Enterobacterales*

Gram positive organisms:

- *Streptococcus* sp.
- Yeasts
- *Staphylococcus* sp.



56

Biochemicals

Conventional use



Modern Use





57

Biochemical Testing

- One of the most important methods for phenotypic microbial identification
- Microbes utilize (or metabolize) biochemicals in unique ways, allowing scientists to accurately identify them
- Examples include:

Carbohydrates	Methyl red
Citric acid	Hydrogen sulfide
Indole	Urea



58

API 20E Strip Biochemicals

- **ONPG**: Detects the enzyme β -galactosidase
- **ADH**: Decarboxylation of the amino acid (aa) arginine by arginine dihydrolase
- **LDC**: Decarboxylation of the aa lysine by lysine decarboxylase
- **ODC**: Decarboxylation of the aa ornithine by ornithine decarboxylase
- **CIT**: Use of citrate as a sole carbon source



59

API 20E Strip Biochemicals

- **H₂S**: Hydrogen sulfide production
- **URE**: Urease enzyme test
- **TDA**: (Tryptophan deaminase): Detection of the tryptophan deaminase enzyme
IND: Indole test: production of indole from tryptophan by the enzyme tryptophanase
- **VP**: The Voges-Proskauer test for the detection of acetoin (acetylmethylcarbinol) produced by fermentation of glucose by bacteria using the butylene glycol pathway



60

API 20E Strip Biochemicals

- **GEL**: Production test of the enzyme gelatinase which liquefies gelatin
- **GLU**: Glucose fermentation (hexose sugar)
- **MAN**: Fermentation of mannose (hexose sugar)
- **INO**: Fermentation of inositol (cyclic polyalcohol)
- **SOR**: Fermentation of sorbitol (sugar alcohol)



API 20E Strip Biochemicals

- **RHA:** Fermentation of Rhamnose (methyl pentose sugar)
- **SAC:** Fermentation of sucrose (disaccharide)
- **MEL:** Fermentation of melibiose (disaccharide)
- **AMY:** Fermentation of amygdalin (glycoside)
- **ARA:** Fermentation of arabinose (pentose sugar)

61



Safety

- Universal precautions: handle all materials as potentially infectious
- Don appropriate PPE at all times – gloves, lab coat, safety glasses
- Utilize aseptic technique
- Inspect packaging before use to ensure its integrity
- Dispose of any materials that come into contact with the specimen or kit components in biohazardous waste
- No food, drink or gum chewing
- Wash hands thoroughly with soap and water after use

62



Materials

- Gloves
- Pipettes
- Ampule protector
- Ampule rack
- Incubation box (tray and lid)
- Deionized water squeeze bottles
- Clip seal
- Package insert
- Scissors
- McFarland demonstration tubes



63



Reagents

- API Test Strip
- API NaCl 0.85% medium
- Mineral oil
- Distilled water



64



Procedural Steps

- Don appropriate PPE
- Prepare the incubation box and distribute 5ml of distilled water into the honey combed wells of the tray to create a humid atmosphere
- Record the identification information on the box
- Place the gallery in the incubation box under sterile conditions



65



Procedural Steps

- Open the ampule of NaCl 0.85%
- Using a swab, remove a single well-isolated colony (18-24 hours old) from the plate
- Carefully emulsify to achieve a homogenous bacterial solution that is 0.5 McFarland
- Immediately use the suspension after preparation (refer to next slide)
- Incubate the tray at $36^{\circ}\text{C} \pm 2^{\circ}\text{C}$ for 18-24 hours




66




Procedural Steps


Using a sterile pipette, fill the compartments with the bacterial suspension as follows:




67



Division of
**Consolidated
Laboratory Services**







DEPARTMENT OF
GENERAL SERVICES



Bacteria and Virus Station 5: API Interpretation

Presenters:
January 23, 2024

68



69

Learning Objectives

- Review the API set-up, briefly
- Observe the final steps of the API strip prior to interpretation
- Analyze the reaction of each API well using the reading table
- Interpret the bacterial identification of the strip using the numerical profile



70

Safety

- Universal precautions: handle all materials as potentially infectious
- Don appropriate PPE at all times – gloves, lab coat, safety glasses
- Utilize aseptic technique
- Inspect packaging before use to ensure its integrity
- Dispose of any material that comes into contact with the specimen or kit components in biohazardous waste
- **No food, drink or gum chewing** at stations
- Wash hands thoroughly after use



Materials

- Gloves
- Pipettes
- Result sheet
- Package insert
- Reading table
- Computer
- APIWEB database (login required)



71



Reagents

- Prepared API strips
- API reagent kit:
 - TDA
 - James
 - VP1 and VP2



72



Procedural Steps – Reading the Strip

- Don appropriate PPE – Gloves, safety glasses
- Record all the reactions on the result sheet
 - Except those that require the addition of more reagents (use the reading table)
- TDA test: Add 1 drop of TDA reagent
 - A reddish brown indicates a positive reaction



73



Procedural steps – Reading the Strip

- VP test: Add 1 drop of VP1 and VP2 reagents. Wait 10 minutes. A pink or red color indicates a positive reaction. Slight pink after 10 minutes should be considered negative
- **Perform the IND test last:** Add 1 drop of JAMES reagent. A pink color developed in the whole cupule indicates a positive reaction



74



Procedural Steps - Reading Table

READING TABLE					
TESTS	ACTIVE INGREDIENTS	QTY (mg/cup)	REACTIONS/ENZYMES	RESULTS	
				NEGATIVE	POSITIVE
ONPG	2-nitrophenyl-β-D-galactopyranoside	0.223	β-galactosidase (Ortho NitroPhenyl-β-D-Galactopyranosidase)	colorless	yellow (1)
ADH	L-arginine	1.9	Arginine DiHydrolase	yellow	red / orange (2)
LDC	L-lysine	1.9	Lysine DeCarboxylase	yellow	red / orange (2)
ODC	L-ornithine	1.9	Ornithine DeCarboxylase	yellow	red / orange (2)
[CIT]	trisodium citrate	0.756	CITrate utilization	pale green / yellow	blue-green / blue (3)
H ₂ S	sodium thiosulfate	0.075	H ₂ S production	colorless / greyish	black deposit / thin line
URE	urea	0.76	UREase	yellow	red / orange (2)
TDA	L-tryptophane	0.38	Tryptophane DeAminase	yellow	TDA / immediate reddish brown
IND	L-tryptophane	0.19	INDole production	JAMES / immediate colorless pale green / yellow	pink
[VP]	sodium pyruvate	1.9	acetoin production (Voges Proskauer)	colorless	VP 1 + VP 2 / 10 min pink / red (5)
[GEL]	Gelatin (bovine origin)	0.6	GELatinase	no diffusion	diffusion of black pigment
GLU	D-glucose	1.9	fermentation / oxidation (GLUcose) (4)	blue / blue-green	yellow / greyish yellow
MAN	D-mannitol	1.9	fermentation / oxidation (MANnitol) (4)	blue / blue-green	yellow
INO	inositol	1.9	fermentation / oxidation (INOsitol) (4)	blue / blue-green	yellow
SOR	D-sorbitol	1.9	fermentation / oxidation (SORbitol) (4)	blue / blue-green	yellow
RHA	L-rhamnose	1.9	fermentation / oxidation (RHAmmose) (4)	blue / blue-green	yellow
SAC	D-sucrose	1.9	fermentation / oxidation (SACcharose) (4)	blue / blue-green	yellow
MEL	D-melibiose	1.9	fermentation / oxidation (MELibiose) (4)	blue / blue-green	yellow
AMY	amygdalin	0.57	fermentation / oxidation (AMYgdalin) (4)	blue / blue-green	yellow
ARA	L-arabinose	1.9	fermentation / oxidation (ARAbinose) (4)	blue / blue-green	yellow
OX	(see oxidase test package insert)		cytochrome-Oxidase	(see oxidase test package insert)	

75



Procedural Steps - Interpretation

- Identification is obtained with a numerical profile
- Look at the results sheet and visualize that the tests are separated into groups of 3 and a value of 1, 2 or 4 is indicated for each.
- Add together the values corresponding to a positive reaction within each group to obtain a 7-digit profile number for the 20 tests
- The oxidase reaction constitutes the 21st test and has a value of 4 if positive

76



Procedural Steps - Interpretation

- If using the database: look up the numerical profile in the list of profiles
- If using the identification software: Enter the 7-digit numerical profile manually into the computer